

#6: Data Acquisition at CBU MEG Lab

Yury Shtyrov



*MRC Cognition and Brain
Sciences Unit*

The main points

- Apologies
- Experimental design: MEG \neq fMRI
- Practicalities: the BB is watching you
- Subject preparation
- Data acquisition *per se*

Triggers

- fMRI: stimuli presented, while scan times are logged, OR, stimuli triggered by the scan
- MEG/EEG: the opposite


MEG Triggers

- Data are recorded *continuously* (1kHz sampling rate = a snapshot of 306 channels every 1 ms)
- Triggers must be presented with you stimuli
- Triggers are recorded with the MEG data
- Response triggers are also fed back to the recorder
- Activity is later assessed in relation to the pre-recorded triggers



Recording data and triggers

activity of the brain
recorded on outside

data processing after (during) the recording



S1 S2 S3 S1 S5 S3 S2 S1



Event-related potentials (ERP)/event-related fields (ERFs)

Timing of stimulation/responses
locked to MEG/EEG signal

‘Trigger-free’ MEG

- Spontaneous data
- E.g. record for 2 min with eyes closed, and for two minutes with eyes open, look at the oscillations in different frequency bands

Exp design

- Think of your stimuli/triggers beforehand
- Better have over- than under-specified stimulation protocol: you can always collapse across a few stim types, but the opposite is more difficult:

1 1 1 1 1 1 1 1 1 2 1 1 2 1 1 1 1 2 1 1 1
1 1 1 1 1 1 1 1 1 2 3 1 2 1 1 1 1 2 3 1 1

Exp design

- Not limited by TR, etc.
- No noise!
- Flexible SOA, length, etc.
- >50 trials for averaging for ER designs (>100 desirable)
- >100 trials for oscillations
- Overall recording time: <2 hrs (ideally ~1h)
 - the subject has to stay still!


Exp design: conditions

- Layers of subtractions... well... maybe not in MEG
- Orthogonal conditions: thinking the ANOVA way

Exp design: conditions

- Orthogonal conditions, an example
- 2x2 design:
 - words attended -words unattended
 - nonwords attended -nonwords unattended
- Data extracted in numerical form from each condition for all individual subject
- ANOVA with two factors: attention X lexicality (plus any additional – e.g. activation locus)

Exp design: conditions



- Data extracted in numerical form from each condition for all individual subject
- ANOVA with two factors: attention X lexicality (plus any additional – e.g. activation locus)
- SPM5 talk by Rik on Jan 29
- dSPM, etc.

Exp design: stimulus presentation

- Auditory (tubes), visual (projection), somatosensory (electric) readily-available
- Software: same stim PC as in the fMRI scanner (E-Prime etc), anything else connectable
- Button-press responses: two four-digit button boxes available

Exp preparation: further considerations

- Baseline: quiet interval before/after stim usually desired
- Subject position: supine vs upright

Wanna do fancy source localisation stuff?

- Get subject with good structural MR images
- Consider recording EEG+MEG

Practical bits

- Get your ethics sorted
- Get your IMC approval
- Follow SOPs
- Obey your MEG operator
- Get the right volunteers

Practical bits

Health and safety issues:

- Liquid Helium
- Working with volunteers

Practical bits

- Liquid Helium: $\sim -270\text{C}$
- Evaporation
- Frost bites
- Asphyxia

Practical bits

Working with volunteers:

- Same basic rules as anywhere, but –
- You will attach things to volunteer skin (coils, electrodes)
- Screen, inform
- Wear gloves, be careful

Practical bits

MEG access:

- Authorised/certified staff
- Only qualified operators can run acquisition
- Other researchers can only be in the MEG lab if accompanied by a qualified operator
- Full de-metallisation on entry

Practical bits

Entrance to the shielded room:

- No metal objects: mobile phones, keys, coins, belts, etc.
- No shoes
- Subjects to be screened

Data acquisition: subject preparation

- Prepare your forms, coils, electrodes, etc. in advance
- De-metallise yourself and your volunteer
- Leave your phone in your office or turn it off
- Brief the subject about the experiment, do the consent, etc.

Data acquisition: subject preparation

- EEG electrodes (1-128, cap)
- EOG (2+2)
- Ground
- HPI coils

Data acquisition: preparation

- Empty room measurement
- ‘Unwired’ subject measurement (esp. for complex set-ups)


Data acquisition: subject preparation

- HPI coils:

4 HPI coils should be placed to maximize distance between them: two on each side of the forehead (right below the hairline) and one behind and above each ear, ensuring that, once the subject is in the machine, all coils will be inside the dewar and not below it. The coils require no paste, and should be affixed using both double-sided sticky rings and tape.

Data acquisition: subject preparation

- Enter subject details in the database (no name, use volunteer ID for name/surname)
- Choose project settings
- Digitise HPI (head position indicator) coils and landmarks
- Digitise additional point/electrodes



Acquisition; control

File On-line Tools Help

name megacq ... 200 Good work, megacq
 wkup ... 200 Connections open

Settings

Change Project: not specified!

Change Subject: Brett Matthew

Change Acquisition: 0.10 ... 330 Hz @ 1000.0 Hz, 321 channels, MaxShield Off

Change External stimulus generation.

Change On-line averaging: not active.

Change HPI: not digitized!

Change Gantry position: Supine

Acquisition controls

GO! Stimulate Average Record raw cHPI

Stopwatch

00:00 Start Stop

No news is good news.
 No raw data recorded.
 Averager not active
 Nothing rejected
 No news from HPI

This is so wonderful

Available subjects :

Subject type : Volunteers

Accessible to group : None

< new >
 Brett Matthew
 Cbu060069 Cbu060069
 Cbu060575 Cbu060575
 Cbu060727 Cbu060727

Subject name :

Last :

First :

Middle :

Date of birth :

day month year

Subject dimensions :

height weight

Other information :

Sex : Male ID # : 0012

Handedness : Right HIS ID :

Comments :

Ok Cancel Help

Coordinate frame alignment

Align frame

Nasion

Status?

x = ?????.? mm
y = ?????.? mm
z = ?????.? mm

LPA

x = ?????.? mm
y = ?????.? mm
z = ?????.? mm



RPA

x = ?????.? mm
y = ?????.? mm
z = ?????.? mm

HPI coils

Coil # 1

x = ?????.? mm
y = ?????.? mm
z = ?????.? mm

Coil # 2

x = ?????.? mm
y = ?????.? mm
z = ?????.? mm

Coil # 3


x = ?????.? mm
y = ?????.? mm
z = ?????.? mm

Coil # 4


x = ?????.? mm
y = ?????.? mm
z = ?????.? mm

EEG not in use

Digitized data



Additional data



No distance computed

OK

Cancel

Help


Digitization

- Have the subject sit in the wooden chair
- Move all objects with large metal components away from the chair and polhemus
- Turn on the polhemus
- Position the goggles on the subject's head, being careful not to cover up any electrodes or HPIs. It is critical that the goggles are tight and do not move once digitization has begun.
- At the top of the screen, hit the “Coordinate frame alignment” button; you are now ready to begin digitization
- digitize the three cardinal landmarks by pointing to each and simultaneously clicking the button on the stylus (electronic pen)
- Check that the two ear points are within a maximum 5 (3?)mm
- Digitise coils
- EEG electrodes, additional points, head shape

Digitization



- Make sure you know the locations!
- Essential for coregistering with MRI images





Data recording

- Done by qualified operators with researcher present all the time
- Recording settings are researcher's responsibility

Acquisition: control

File On-line Tools Help

name megacq ... 200 Good work, megacq
wkup ... 200 Connections open

Settings

Change Project: not specified!

Change Subject: Brett Matthew

Change Acquisition: 0.10 ... 330 Hz @ 1000.0 Hz, 321 channels, MaxShield Off

Change External stimulus generation.

Change On-line averaging: not active.

Change HPI: not digitized!

Change Gantry position: Supine

Acquisition controls

GO! Stimulate Average Record raw cHPI

Stopwatch

00:00 Start Stop

No news is good news.
No raw data recorded.
Averager not active
Nothing rejected
No news from HPI

Select a project

Available projects :

- < new >
- camsomstim
- camtest
- mmneeg
- motortest

Identifier :

ID # : 0010

Name :

Aims of this project :

Responsible persons :

Comments :

OK Cancel Help

How to acquire?

Channel selection :

MEG 0113	MEG 0112	MEG 0111	MEG 0122	MEG 0123	MEG 0121	MEG 0132	MEG 0133	MEG 0131	MEG 0143	MEG 0142	MEG 0141	MEG 0213	MEG 0212	MEG
MEG 0313	MEG 0312	MEG 0311	MEG 0322	MEG 0323	MEG 0321	MEG 0333	MEG 0332	MEG 0331	MEG 0343	MEG 0342	MEG 0341	MEG 0413	MEG 0412	MEG
MEG 0513	MEG 0512	MEG 0511	MEG 0523	MEG 0522	MEG 0521	MEG 0532	MEG 0533	MEG 0531	MEG 0542	MEG 0543	MEG 0541	MEG 0613	MEG 0612	MEG
MEG 0713	MEG 0712	MEG 0711	MEG 0723	MEG 0722	MEG 0721	MEG 0733	MEG 0732	MEG 0731	MEG 0743	MEG 0742	MEG 0741	MEG 0813	MEG 0812	MEG
MEG 0932	MEG 0933	MEG 0931	MEG 0942	MEG 0943	MEG 0941	MEG 1013	MEG 1012	MEG 1011	MEG 1023	MEG 1022	MEG 1021	MEG 1032	MEG 1033	MEG
MEG 1133	MEG 1132	MEG 1131	MEG 1142	MEG 1143	MEG 1141	MEG 1213	MEG 1212	MEG 1211	MEG 1223	MEG 1222	MEG 1221	MEG 1232	MEG 1233	MEG
MEG 1333	MEG 1332	MEG 1331	MEG 1342	MEG 1343	MEG 1341	MEG 1412	MEG 1413	MEG 1411	MEG 1423	MEG 1422	MEG 1421	MEG 1433	MEG 1432	MEG
MEG 1533	MEG 1532	MEG 1531	MEG 1543	MEG 1542	MEG 1541	MEG 1613	MEG 1612	MEG 1611	MEG 1622	MEG 1623	MEG 1621	MEG 1632	MEG 1633	MEG
MEG 1732	MEG 1733	MEG 1731	MEG 1743	MEG 1742	MEG 1741	MEG 1813	MEG 1812	MEG 1811	MEG 1822	MEG 1823	MEG 1821	MEG 1832	MEG 1833	MEG
MEG 1932	MEG 1933	MEG 1931	MEG 1943	MEG 1942	MEG 1941	MEG 2013	MEG 2012	MEG 2011	MEG 2023	MEG 2022	MEG 2021	MEG 2032	MEG 2033	MEG
MEG 2133	MEG 2132	MEG 2131	MEG 2143	MEG 2142	MEG 2141	MEG 2212	MEG 2213	MEG 2211	MEG 2223	MEG 2222	MEG 2221	MEG 2233	MEG 2232	MEG
MEG 2332	MEG 2333	MEG 2331	MEG 2343	MEG 2342	MEG 2341	MEG 2412	MEG 2413	MEG 2411	MEG 2423	MEG 2422	MEG 2421	MEG 2433	MEG 2432	MEG
MEG 2533	MEG 2532	MEG 2531	MEG 2543	MEG 2542	MEG 2541	MEG 2612	MEG 2613	MEG 2611	MEG 2623	MEG 2622	MEG 2621	MEG 2633	MEG 2632	MEG
STI 7	STI 8	STI 9	STI 10	STI 11	STI 12	STI 13	STI 14	STI 15	STI 16	EEG 1	EEG 2	EEG 3	EEG 4	EEG
EEG 15	EEG 16	EEG 17	EEG 18	EEG 19	EEG 20	EEG 21	EEG 22	EEG 23	EEG 24	EEG 25	EEG 26	EEG 27	EEG 28	EEG
EEG 39	EEG 40	EEG 41	EEG 42	EEG 43	EEG 44	EEG 45	EEG 46	EEG 47	EEG 48	EEG 49	EEG 50	EEG 51	EEG 52	EEG
EEG 63	EEG 64	EEG 65	EEG 66	EEG 67	EEG 68	EEG 69	EEG 70	EEG 71	EEG 72	EEG 73	EEG 74	EEG 75	EEG 76	EEG
EEG 87	EEG 88	EEG 89	EEG 90	EEG 91	EEG 92	EEG 93	EEG 94	EEG 95	EEG 96	EEG 97	EEG 98	EEG 99	EEG 100	EEG
EEG 111	EEG 112	EEG 113	EEG 114	EEG 115	EEG 116	EEG 117	EEG 118	EEG 119	EEG 120	EEG 121	EEG 122	EEG 123	EEG 124	EEG
MISC 7	MISC 8	MISC 9	MISC 10	MISC 11	MISC 12	MISC 13	MISC 14	MISC 15	MISC 16	STI 101	STI 102	STI 201	STI 301	MISC
MISC 205	MISC 206	MISC 301	MISC 302	MISC 303	MISC 304	MISC 305	MISC 306							

Number of channels = 321

OK

Cancel

Help

Event characteristics

Event: 1 Channel: STI101 New: 1/63 Old: 0/63 Delay: 0 ms

Event: 2 Channel: STI101 New: 2/63 Old: 0/63 Delay: 0 ms

Event: 3 Channel: STI101 New: 4/63 Old: 0/63 Delay: 0 ms

Event: 4 Channel: STI101 New: 8/63 Old: 0/63 Delay: 0 ms

Event: 5 Channel: STI101 New: 16/63 Old: 0/63 Delay: 0 ms

Averaging categories

01 On-line display t = -50.0.. 300.0 ms Event: 1 No condition Sub: 0 No limit for # of averages

02 On-line display t = -50.0.. 300.0 ms Event: 1 No condition Sub: 0 No limit for # of averages

03 On-line display t = -50.0.. 300.0 ms Event: 1 No condition Sub: 0 No limit for # of averages

Artefact rejection settings

	MEG grad	MEG mag	EEG	EOG	EMG	ECG
Max :	<input type="text" value="3000"/> fT/cm	<input type="text" value="3000"/> fT	<input type="text" value="1"/> uV	<input type="text" value="150"/> uV	<input type="text" value="1"/> uV	<input type="text" value="1"/> uV
Slope :	<input type="text" value="1"/> fT/cm	<input type="text" value="1"/> fT	<input type="text" value="1"/> uV			
Spike :	<input type="text" value="1"/> fT/cm	<input type="text" value="1"/> fT	<input type="text" value="1"/> uV			
No signal :	<input type="text" value="1"/> fT/cm	<input type="text" value="1"/> fT				
Noisy :	<input type="text" value="1"/> fT/cm	<input type="text" value="1"/> fT				
Ignore (ms) after stimulus :	<input type="text" value="0"/>					

Miscellaneous settings

Display update interval (s) :

Data recording


- Record raw/online averaging only
- Set-up stim categories for online averaging
 - Stim types
 - Epoch length
 - Baseline
 - Subaverages
 - Artifact rejection

Data recording

- Record raw/online averaging only
- Set-up stim categories for online averaging
- Sampling rate (1kHz)
- Bandpass
- Channels (EEG, EOG, stim, misc)
- Maxshield
- Continuous HPI recording
- Save settings – you can load them later (more than 1 possible)


Data recording: all ready

- De-metallise
- Take the subject in
- Connect EEG, EOG, ground, HPI coils
- Ear pieces/button boxes/video etc (check in advance)
- Move the subject up as close to the helmet and possible, make sure they are comfortable
- Intercom/camera
- Lock'em up and start the show!



- **Go!**
- **Instruct the subject to stay still**
- **Measure HPI**
- **Select Average and/or Record Raw**
- **cHPI**

Acquisition: raw data





▲ ▼ || Selection... Rack Scales... Colors... xy Window: 3.0 s

[MEG 01-02] Buffer 18.00 s...

The end

- Monitor channels
- Monitor rejections
- Monitor averages
- Monitor subject
- Give them breaks, talk to them
- Stop – save averages
AND raw data
- Repeat for each block
(15-20 min)





Time integration Calculate slopes No reference cursor

autoscale

visualleftrightupdown.fif

Next week

Lecture 7:

From signals to event-related fields and dipoles: an introduction using Elekta-Neuromag software

15.1.2007

THE END

