


64-channel electrode configuration for MEGIN Triux neo at CBU



- 1 PO9
- 2 FPz
- 3 PO10
- 4 AF7
- 5 AF3
- 6 Afz
- 7 AF4
- 8 AF8
- 9 F7
- 10 F5
- 11 F3
- 12 F1
- 13 Fz
- 14 F2
- 15 F4

16 F6
17 F8
18 FT7
19 FC5
20 FC3
21 FC1
22 FCz
23 FC2
24 FC4
25 FC6
26 FT8
27 T7
28 C5
29 C3
30 C1
31 Cz
32 C2
33 C4
34 C6
35 T8
36 TP9
37 TP7
38 CP5
39 CP3
40 CP1
41 CPz
42 CP2
43 CP4
44 CP6
45 TP8
46 TP10
47 P7
48 P5
49 P3
50 P1
51 Pz
52 P2
53 P4
54 P6
55 P8
56 PO7
57 PO3
58 Poz
59 PO4
60 PO8
61 O1

62 Oz

63 O2

64 Iz